IC1022

MASSAGE CHAIR
FAUTEUIL DE MASSAGE

USER MANUAL

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MASSAGE CHAIR
FAUTEUIL DE MASSAGE

SAVE THE INSTRUCTIONS

Thank you for purchasing this product. Please read this instruction manual carefully before operating this massage chair.

Please pay attention to the important safety warnings and keep this manual for future use. Note: EBI reserves the right to revise the design and description of this product without notice.

CONTENTS

Safety and maintenance	2
Name and function of components	4
Features.....	5
Methods of use	
1. The use of the controls	6
2. Characteristics of automatic massage features	6
3. Auto---timing is 15 minutes	8
4. Manual function	9
5. Mp3 display	11
6. Calf rest and backrest recline and decline function	11
7. Method to dowload mp3 music.....	12
8. Padding and back cushion assembly and operation	14
9. Installation and use of the controller and tea table	14
10. Improper use	15
11. Floor protection	16
Product specifications.....	16
FRANÇAIS	19

SAFETY AND MAINTENANCE

1. IMPORTANT SAFETY WARNING

- Do not allow children to touch moveable part of this product.
- Please use a grounded power supply suited to this product.
- Please pull out the plug after use or before cleaning to avoid injury or damage to the product.
- Please operate this product according to the instructions In this manual.
- Do not use accessories which are not recommended by EBI.
- Do not use this product outdoors.
- Please read this instruction manual carefully before operation.
- Recommended use time is 15 minutes.
- Please do not use the chair if the leather is broken or damaged.
- Please do not use it if the radiator hole is covered or blocked.
- Please do not drop anything into this product.
- Please do not fall asleep while using this product.
- Please do not use the chair while under the influence of drugs or alcohol.
- Please do not use the chair within one hour after having meal.

2. ENVIRONMENT FOR USE

- Please do not use the chair under high temperature or moist environment such as bathroom.
- Please do not use the chair immediately if the environmental temperature changes sharply.
- Please do not use the chair in dusty environments.
- Please do not use the chair in a constrained space or poorly ventilated areas.

3. PEOPLE WHO SHOULD REFRAIN FROM USING THIS PRODUCT

- People suffering from osteoporosis should refrain from using this product.
- People suffering from heart disease or carrying electronic medical equipment such as a heart pacemaker should refrain from using this product.
- People suffering from fever should refrain from using this product.
- Pregnant women should refrain from using this product.
- People with injuries or skin disease should refrain from using this product.
- Children under twelve or people who are not able to operate this product should refrain from using this product.

4. SAFETY

- Check that the voltage is suitable to the specifications of this product.
- Do not use wet hands to pull the plug.
- Do not spill liquids onto this product to avoid electric shock or cause damage to this product.
- Do not damage the wires or change the circuits of this product.
- Do not use a wet cloth to clean the electric parts such as switch and plug.
- In the event of a power failure un-plug the chair
- Stop using this product while it works abnormal, and please, consult local agent immediately.
- Stop use of this product if you feel unwell, and consult a physician.

5. MAINTENANCE


- This product should be maintained by certified technicians.
- DO NOT attempt to disassemble or maintain this chair.
- Disconnect the power cord after use.
- Do not use this product if the socket is loose.
- If this product will be left unused for long periods of time,
- Store this product in dry and dustless environment.
- Do not store this product under high temperature or near fire and avoid
- Long exposure to direct sunshine.
- Clean this product with a dry cloth. Do not use thinner, benzene or alcohol.
- The mechanical components of this product are specially designed and
- Do not require maintenance.
- Do not move the chair using the wheels on wood floor, rough ground or in a narrow space. Two people should lift the chair together to prevent damage.
- Please use this chair intermittently; do not use this product continuously for long periods of time.

6. TROUBLESHOOTING


- It is normal for the motor to make noise while in use.
- If the controller doesn't work normally, please check that the plug and socket are properly connected and the switch is turned on.
- When the rated working time is over, the product will automatically activate the power-off button; If this product works continuously for long, the temperature protector will make the product automatically power-off.
- It can be used again after half an hours rest.

NAME AND FUNCTION OF COMPONENTS

1. Pillow Pad
2. Remote Control
3. Earphone
4. Armrest
5. Cushion
6. Calfrest
7. Footrest
8. Backcushion
9. Arms airbag
10. Jade Heater
11. Mini Tray

12. Backrest Cover
13. Fan Hole
14. Remote Controller Stand
15. Minitray Stand(Hidden)
16. MP3 Player
17. Back Wheel
18. Earphone Line
19. Cable & Plug
20. Earphone Jack
21. Plug Jack
22. Switch
23. Fusebox

FEATURES

- Dispel fatigue.
- Accelerate blood circulation.
- Alleviate muscle fatigue.
- Relax tense muscles.
- Alleviate nerve aches.
- Alleviate muscle aches.
- Massage mechanism is a four-wheel driven, muted design and can move Up & down.
- Shiatsu, kneading, tapping, knocking, kneading & tapping and rolling six emulational massage modes
- Automatically checks the upper body massage points & shoulder position.
- Five characteristic massage functions (relax, health care, physical therapy, fatigue and prescription).
- Upper body automatic massage (neck & shoulder, whole upper body, back & waist, overall, partial and fixed - six choices).
- Under fixed state, the massage mechanism can move up and down to reach an accurate position.
- Lower body airbag massage function: equipped with 33 airbags, three modes and three levels of speed. (Strong, medium and weak).
- Arms airbag massage function: equipped with 36 airbags.
- Vibration massage function: set in the back and buttock position, there are three choices: strong, weak and off.
- Music display and volume adjustment function; previous song, next song choices,
- Music accompaniment function: back massage speed and vibration changes with the high frequency, buttock vibration strength changes with the middle frequency.
- Under the music accompaniment state, the controller in the wire can be used to adjust the massage strength.
- Automatic recline and decline function: it is designed with modes that the calf rest and backrest recline and decline separately, lie down, automatically or manually recover after turning off the chair.
- The massage mechanism automatically recovers when turning off the chair.
- The chair equipped with a heater which uses five Xiu jade as the infrared emission source, the temperature can be adjusted between 40 to 70 degrees celsius.
- It is equipped with a seven color LCD, MP3 display, earphones and mini tea table.

METHODS OF USE

1. THE USE OF THE CONTROLS

- **Start/Stop massage:** press power button to start massage. During the massage period, press the power button and it will stop, the calf rest will return to its original position(move down to the bottom)
- **Neck position adjustment:** press up and down button

2. CHARACTERISTICS OF AUTOMATIC MASSAGE FEATURES

- **Relax:** massage the whole upper body and also start all other functions (lower body airbag and vibration massage) to relax muscle. Set the working time to be 8 minutes. Mode in display window shows (01).

- **Health care:** besides the whole upper body massage function, it can massage respectively shoulders, back and waist. You can also start all other functions to accelerate blood circulation. Set the working time to be 15 minutes. Mode in display window shows (02).

- **Physical therapy:** massage shoulder, back and waist, in the meantime start all other functions to relax aches in the shoulder, back and waist. Set the working time to be 30 minutes. Mode in display window shows (03).

- **Fatigue:** deeply massage whole upper body, shoulder, back and waist. In the meantime start all other functions to relax whole body. Set the working time to be 38 minutes. Mode display window shows (05).

- **Personal:** A total of 16 massage functions. The upper body is divided into 16 massage points. Enjoy a deep massage according to your preferences. Relax specific muscles aches aches.

- **Awake:** automatically start massage function. During the massage period, press (awake) button then all massage functions will stop and the timing display window shows (01), then press again the (awake) button to set the automatic start time. The automatic start time can also be set when the massage chair is turned off. The max timing is (08) (that is 8 hours). After setting the automatic start time, you may lie on the chair to sleep or relax. Half an hour before the automatic start time ends, the chair will begin to automatically remind you that time is up. Timing display window shows the automatic start time. Mode in display window shows (99)

E.g.: press (awake) button continually three times, then the timing display window shows (03) which means the automatic start time is 3 hours, and the mode display window shows (99) which means the awake function is turned on. From the moment you finished setting the time, two hours and a half later the massage function will start automatically.

3. AUTO---TIMING IS 15 MINUTES.

- **Whole upper body program**---shows (F0).

- **Partial program**---shows (F1). Press up and down button to move the mechanism to the required position.

- **Fixed program**---shows (F2) Press up and down button to move the mechanism to the required position.

- **Shoulder program**---shows (F3).

- **Back program**--shows (F4).

- **Waist program**--shows (F5).

4. MANUAL FUNCTION.

- **Timing:** it can only be set under manual massage function; it is designed with three levels: 15 minutes, 20minutes and 30minutes.
- **Manual:** total 6 manual choices: fixed, partial, overall, shoulder, back and waist shows respectively P0, P1, P2, P3, P4 and P5.
- **Massage methods:**
Tapping, knocking, kneading & tapping, rolling,shiatsu and kneading.

- **Width:** adjustment of the width of massage mechanism. It is designed with three levels: wide, middle and narrow.

- **Speed:** adjustment of the massage speed. The speed can't be changed under the shiatsu mode, and under other modes. There are five levels.

- **Back vibration:** back vibration massage strength and on/off button.

- **Buttock vibration:** buttock vibration strength and on/off button.

- **Arms:** arms airbag strength and on/off button.

- **Lower body:** lower body airbag on/off button.

- **Strength:** lower body airbag strength button, three levels: strong, medium and weak.

- **Mode:** lower body airbag modes.

There are three modes: Lower body, Foot & Calf, Buttocks.

- **Handheld heat therapy:** on/off button of the handheld heat therapy. Default temperature is 40 . (there are seven levels available: 40 45 50 55 60 65 70).

- **Handheld heat therapy (+):** increase the temperature of the handheld heater. It will increase by one level upon each press, the highest temperature is 70.

- **Handheld heat therapy (-):** decrease the temperature of the handheld heater. It will decrease by one level upon each press. The lowest temperature is 40.

- **Background color:** the background color changes automatically with the different massage methods.
- **Color adjustment:** manually set the background color. There are seven colors: red, green, blue, purple, orange, white and cyan.

5. MP3 DISPLAY

- **Functions:** turn on, turn off, previous song, next song, synchronization.

- **Music accompaniment function:**
 1. Back massage speed and vibration strength changes with the high frequency.
 2. Buttock vibration strength changes with the medium frequency.

6. CALF REST AND BACKREST RECLINE AND DECLINE FUNCTION.

It is designed with calf rest recline and decline, backrest recline and decline, sit and lie down.


7. METHOD TO DOWNLOAD MP3 MUSIC

- 1) Turn off the power of the massage chair and open the cover of the MP3 box behind the massage chair.
- 2) Gently press the SD card and it will pop out of the MP3 player, then remove it.
- 3) Insert the SD card into the SD card read-write equipment supplied with the chair.
- 4) Insert the SD card read-write equipment into the USB port of your PC.(windows XP, windows 2000, windows NT)
- 5) Download your favorite song in MP3 format onto the SD card (movable disk).
- 6) Remove the USB Mass Storage Device (movable disk) from the PC safely.
- 7) Pull the SD card out of the USB port.
- 8) Pull the SD card out of the SD card read-write equipment.
- 9) Insert the SD card into the SD card port of the MP3 player.
- 10) The refresh of the MP3 music is done, can turn on the power of the massage chair again.
- 11) If the MP3 can not be played, please format the SD card in the PC, download MP3 type music, then insert it into the massage chair as above.

ATTENTION: DO NOT REMOVE OR DISASSEMBLE THE MP3 PLAYER. ONLY THE SD CARD CAN BE REMOVED.

- (1) MP3 box
- (2) SD card
- (3) front cover of the SD card read-write equipment.
- (4) SD card read-write equipment
- (5) SD card
- (6) back cover of the SD card read-write equipment.


8. PADDING AND BACK CUSHION ASSEMBLY AND OPERATION

The removable cushion can reduce the neck & shoulder kneading Massage. You can determine according to your demands whether to use the pillow cushion or not (We recommend using the pillow cushion).


The back cushions connect to the the zipper(1) The pillow cushion connects to the back Cushion with the velcro strip(2)


Warning: DO NOT raise the back cushion when in use. This may cause injury.

9. INSTALLATION AND USE OF THE CONTROLLER AND TEA TABLE

Install the controller shelf on the massage chair(1). Insert the controller in the controller shelf(2). The control can rotate vertically(3). The control can rotate horizontally(4). Install the tea table tray on the massage chair (5) Put the tea table tray on the shelf (6).

10. IMPROPER USE:

- When the arm airbag is connected to the armrest, DO NOT lubricate the spindle.
- DO NOT stab the airbag with needles or other sharp objects.
- DO NOT wobble the arm airbags. This could damage the air injection tubes.
- When the arm airbag spindle is inserted into the hole half length and the chair is turned on, the induction equipment in the hole receives a signal and will start feeding the air; when the spindle is inserted into the hole completely, there will be no air coming out of the joint and the airbag will start working.
- When the airbag is taken off but not completely, it is normal that air comes out of the joint.
- If you remove the airbag assembly, verify that it moves forward and backwards freely.
- In order to avoid colliding with the airbag assembly, the controller should be tilted forward by 20 degrees as shown in pictures (1) and (2).


11. FLOOR PROTECTION

Prolonged use of this massage chair can damage wood floors. In order to prevent permanent damage, place a carpet or mat under the chair.

CAUTION: DO NOT MOVE THE CHAIR USING THE WHEELS ON WOOD FLOOR, ROUGH GROUND OR IN A NARROW SPACE. TWO PEOPLE SHOULD LIFT THE CHAIR TOGETHER TO PREVENT DAMAGE.

PRODUCT SPECIFICATIONS

Input Voltage: AC110-120V 60Hz, AC220-240V 50Hz

Input power: 220W

**WARRANTY CARD**

ART. NAME		DEALER(SERVICE CENTER)
MODEL		
CUSTOMER		
ADDRESS		
TEL		
DATE OF PURCHASE		EXPIRY : ONE YEAR FROM DATE OF PURCHASE

If the chair malfunctions within the warranty period (one year) and under the following conditions, this product will be repaired for free:

1. there are broken parts or malfunction when purchasing this machine.
2. damage that occurs when it is being used under the guidance of the operation instructions (but the manufacturer is not responsible for damage resulting from misuse, operation without the guidance of the instructions, or natural disasters)
3. the damage of the nonconsumptive parts.

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Pictures are for reference only, colors and models may vary without notice. Made in China. iComfort™ is a trademark of EBI Inc. Photographies pour référence seulement. Les couleurs et modèles sont sujet à varier sans préavis. Fabriqué en Chine. iComfort™ est une marque de commerce de EBI inc.